

README.TXT

Below find a list by product of action done by CSSecurityTool:

Windows Defender: List of folders and executables add in exclusions

Windows Firewall: List of allowed ports

Windows Security: List of Carestream folders with needed rights for the current user

CSImaging7

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Carestream\CS3DSuite-light-uninst.exe

C:\Program Files\Carestream\3D Visualization Application\3DImagingSoftware.exe

C:\Program Files\Carestream\3D Visualization Application\3DViewer-uninst.exe

C:\Program Files\Carestream\3D Visualization Application\Setup3DViewer.exe

C:\Program Files\Carestream\3D Visualization Application\ViewData.exe

C:\Program Files\Carestream\Film Composer\FilmComposer.exe

C:\Program Files\Carestream\Film Composer\storescu.exe

C:\Program Files\Carestream\Film Composer\UninstallFilmComposer.exe

C:\Program Files\Carestream\Quick Burn\QuickBurn.exe
C:\Program Files\Carestream\Quick Burn\Uninstall-QuickBurn-Carestream.exe
C:\Program Files\Carestream\Report Editor\ReportEditor-uninst.exe
C:\Program Files\Carestream\Report Editor\ReportEditor.exe
C:\Program Files\Carestream\Volume Converter\VolumeConverter-uninst.exe
C:\Program Files\Carestream\Volume Converter\VolumeConverter.exe
C:\Program Files (x86)\Carestream\CSI Suite uninst.exe
C:\Program Files (x86)\Carestream\CSDML\AcqTwain.exe
C:\Program Files (x86)\Carestream\CSDML\codebase_cli.exe
C:\Program Files (x86)\Carestream\CSDML\CSDMLite.exe
C:\Program Files (x86)\Carestream\CSDML\CSDMLwdog.exe
C:\Program Files (x86)\Carestream\CSDML\DIStoCSDMLite.exe
C:\Program Files (x86)\Carestream\CSDML\export_ps.exe
C:\Program Files (x86)\Carestream\CSDML\fttfx.exe
C:\Program Files (x86)\Carestream\CSDML\shutdown_tool.exe
C:\Program Files (x86)\Carestream\CSDML\uninstallCSDML.exe
C:\Program Files (x86)\Carestream\CSDML\DcmToDis\DcmToDis.exe
C:\Program Files (x86)\Carestream\CSDML\DcmToDis\DcmToDisUnsafe.exe
C:\Program Files (x86)\Carestream\CSImaging\CSI-uninst.exe
C:\Program Files (x86)\Carestream\CSImaging\FMSEditor.exe
C:\Program Files (x86)\Carestream\CSImaging\LicenseTWDlg.exe
C:\Program Files (x86)\Carestream\CSImaging\RJViewer.exe
C:\Program Files (x86)\Carestream\CSImaging\tc.exe
C:\Program Files (x86)\Carestream\CSImaging\TConvert.exe
C:\Program Files (x86)\Carestream\CSImaging\TrophyPreview.exe
C:\Program Files (x86)\Carestream\CSImaging\ts.exe
C:\Program Files (x86)\Carestream\CSImaging\TW.exe
C:\Program Files (x86)\Carestream\Film Composer\FilmComposer.exe
C:\Program Files (x86)\Carestream\Film Composer\storescu.exe
C:\Program Files (x86)\Carestream\Film Composer\UninstallFilmComposer.exe
C:\Program Files (x86)\Carestream\Monitor\monitor.exe

C:\Program Files (x86)\Carestream\Monitor\UninstallMonitor.exe

C:\Program Files (x86)\Carestream\Patient Browser\DicomConvertToFp.exe

C:\Program Files (x86)\Carestream\Patient Browser\DICOMDIRTool.exe

C:\Program Files (x86)\Carestream\Patient Browser\Patient.exe

C:\Program Files (x86)\Carestream\Patient Browser\uilog.exe

C:\Program Files (x86)\Carestream\Patient Browser\uninst.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\CSAcqDmn.exe

C:\Program Files (x86)\Common Files\Trophy\CSMeshRenderer\uninst.exe

C:\Program Files (x86)\Common Files\Trophy\Installs\Light_SetupRA_SDK.exe

C:\Program Files (x86)\Common Files\Trophy\Installs\Light_Trophy_Acquisition_SDK_Install.exe

C:\Program Files (x86)\Common
Files\Trophy\Installs\Light_Trophy_Processings2D_SDK_with_Tools.exe

C:\Program Files (x86)\Common Files\Trophy\IPET\bin\ipet.exe

C:\Program Files (x86)\Common Files\Trophy\IPET\bin\uninst.exe

C:\Program Files (x86)\Common Files\Trophy\Licensing\activate.exe

C:\Program Files (x86)\Common Files\Trophy\Licensing\uninst.exe

C:\Program Files (x86)\Common Files\Trophy\Mesh Converter\CSMeshConverter.exe

C:\Program Files (x86)\Common Files\Trophy\Mesh Converter\UninstallCSMeshConverter.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomConvertToFp.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomImg.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomPack.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomPdf.exe

C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DEclipse\EclipseUninstaller.exe

C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DPresets\Factory\FactoryPresetUninstaller.exe

C:\Program Files (x86)\Common Files\Trophy\RA_SDK\Appli DemoRA_SDK.exe

C:\Program Files (x86)\Common Files\Trophy\RA_SDK\Light_RA_SDKUninstaller.exe

C:\Program Files (x86)\Common
 Files\Trophy\SDK\Processings2D\bin\win32\release\Processings2DSampleAdvanced.exe
 C:\Program Files (x86)\Common
 Files\Trophy\SDK\Processings2D\bin\win32\release\Processings2DSampleBasic.exe
 C:\Program Files (x86)\Common Files\Carestream\Processing\EclipseDBUninstall.exe
 C:\Program Files (x86)\Common
 Files\Carestream\Processing\Eclipse\ImPrEclipse\CreatePrImage.exe
 C:\Program Files (x86)\Common Files\Carestream\Processing\Eclipse\ImPrEclipse\EclipseProc.exe
 C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe
 C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe
 C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe
 C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe

Windows Firewall

9999	TCP	In/out
10001-10100	TCP	In/out
50000	TCP	In/out
52540	TCP	In/out
52541	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream
 C:\Program Files\Common Files\Trophy
 C:\Program Files (x86)\Carestream
 C:\Program Files (x86)\Common Files\Carestream
 C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CSImaging8_v2

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files (x86)\Carestream\CSI Suite uninst.exe

C:\Program Files (x86)\Carestream\CSDML\AcqTwain.exe

C:\Program Files (x86)\Carestream\CSDML\command_tool.exe

C:\Program Files (x86)\Carestream\CSDML\CSDMLewdog.exe

C:\Program Files (x86)\Carestream\CSDML\CSDMLite.exe
C:\Program Files (x86)\Carestream\CSDML\export_ps.exe
C:\Program Files (x86)\Carestream\CSDML\fttfx.exe
C:\Program Files (x86)\Carestream\CSDML\set_configuration.exe
C:\Program Files (x86)\Carestream\CSDML\shutdown_tool.exe
C:\Program Files (x86)\Carestream\CSDML\sqlite3.exe
C:\Program Files (x86)\Carestream\CSDML\uninstallCSDML.exe
C:\Program Files (x86)\Carestream\CSDML\IngestCodebase\codebase_cli.exe
C:\Program Files (x86)\Carestream\CSDML\IngestCodebase\ingestcodebase.exe
C:\Program Files (x86)\Carestream\CSImaging8\AcqTwain.exe
C:\Program Files (x86)\Carestream\CSImaging8\AnnotationRenderer.exe
C:\Program Files (x86)\Carestream\CSImaging8\csds.exe
C:\Program Files (x86)\Carestream\CSImaging8\csds_config.exe
C:\Program Files (x86)\Carestream\CSImaging8\QuickBurn.exe
C:\Program Files (x86)\Carestream\CSImaging8\ReconstructionHandler.exe
C:\Program Files (x86)\Carestream\CSImaging8\TConvert.exe
C:\Program Files (x86)\Carestream\CSImaging8\TW.exe
C:\Program Files (x86)\Carestream\CSImaging8\TWV.exe
C:\Program Files (x86)\Carestream\CSImaging8\uilog.exe
C:\Program Files (x86)\Carestream\CSImaging8\uninst.exe
C:\Program Files (x86)\Carestream\CSImaging8\UninstallTConvert.exe
C:\Program Files (x86)\Carestream\CSImaging8\FMSEditor\FMSEditor.exe
C:\Program Files (x86)\Carestream\CSImaging8\RJViewer\RJViewer.exe
C:\Program Files (x86)\Carestream\Film Composer\FilmComposer.exe
C:\Program Files (x86)\Carestream\Film Composer\storescu.exe
C:\Program Files (x86)\Carestream\Film Composer\UninstallFilmComposer.exe
C:\Program Files (x86)\Carestream\Monitor\monitor.exe
C:\Program Files (x86)\Carestream\Monitor\UninstallMonitor.exe
C:\Program Files (x86)\Common Files\Carestream\Processing\EclipseDBUninstall.exe
C:\Program Files (x86)\Common
Files\Carestream\Processing\Eclipse\ImPrEclipse\CreatePrImage.exe

C:\Program Files (x86)\Common Files\Carestream\Processing\Eclipse\ImPrEclipse\EclipseProc.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files (x86)\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files (x86)\Common Files\Trophy\ExtraoralFms\extraoral_fms.exe

C:\Program Files (x86)\Common Files\Trophy\ExtraoralFms\uninst.exe

C:\Program Files (x86)\Common Files\Trophy\IPET\bin\ipet.exe

C:\Program Files (x86)\Common Files\Trophy\IPET\bin\uninst.exe

C:\Program Files (x86)\Common Files\Trophy\Licensing\activate.exe

C:\Program Files (x86)\Common Files\Trophy\Licensing\Uninstall activate.exe

C:\Program Files (x86)\Common Files\Trophy\Licensing\resources\elevate.exe

C:\Program Files (x86)\Common Files\Trophy\Mesh Converter\CSMeshConverter.exe

C:\Program Files (x86)\Common Files\Trophy\Mesh Converter\UninstallCSMeshConverter.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomConvertToFp.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomImg.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomPack.exe

C:\Program Files (x86)\Common Files\Trophy\Processing\DicomPdf.exe

C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DEclipse\EclipseUninstaller.exe

C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DEquinox\EquinoxUninstaller.exe

C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DPresets\Factory\FactoryPresetUninstaller.exe

C:\Program Files (x86)\Common Files\Trophy\RA_SDK\Appli DemoRA_SDK.exe

C:\Program Files (x86)\Common Files\Trophy\RA_SDK\Light_RA_SDKUninstaller.exe

C:\Program Files (x86)\Common
Files\Trophy\SDK\Processings2D\bin\win32\release\Processings2DSampleAdvanced.exe

C:\Program Files (x86)\Common
Files\Trophy\SDK\Processings2D\bin\win32\release\Processings2DSampleBasic.exe

C:\Program Files (x86)\Common Files\Trophy\Trace\Trace.exe

C:\Program Files (x86)\Common Files\Trophy\Trace\UninstallTrace.exe

C:\Program Files\Carestream\Volume Converter\VolumeConverter-uninst.exe

C:\Program Files\Carestream\Volume Converter\VolumeConverter.exe
C:\Program Files\Carestream\CS3DSuite-light-uninst.exe
C:\Program Files\Carestream\3D Visualization Application\3DImagingSoftware.exe
C:\Program Files\Carestream\3D Visualization Application\3DViewer-uninst.exe
C:\Program Files\Carestream\3D Visualization Application\Setup3DViewer.exe
C:\Program Files\Carestream\3D Visualization Application\ViewData.exe
C:\Program Files\Carestream\Film Composer\FilmComposer.exe
C:\Program Files\Carestream\Film Composer\storescu.exe
C:\Program Files\Carestream\Film Composer\UninstallFilmComposer.exe
C:\Program Files\Carestream\Quick Burn\QuickBurn.exe
C:\Program Files\Carestream\Quick Burn\Uninstall-QuickBurn-Carestream.exe
C:\Program Files\Carestream\Report Editor\ReportEditor-uninst.exe
C:\Program Files\Carestream\Report Editor\ReportEditor.exe

Windows Firewall

9999	TCP	In/out
10001-10100	TCP	In/out
50000	TCP	In/out
52540	TCP	In/out
52541	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream
C:\Program Files\Common Files\Trophy
C:\Program Files (x86)\Carestream
C:\Program Files (x86)\Common Files\Carestream
C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CSImaging8_v3

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Carestream\CS3DSuite-light-uninst.exe

C:\Program Files\Carestream\3D Visualization Application\3DImagingSoftware.exe

C:\Program Files\Carestream\3D Visualization Application\3DViewer-uninst.exe

C:\Program Files\Carestream\3D Visualization Application\Setup3DViewer.exe

C:\Program Files\Carestream\3D Visualization Application\ViewData.exe

C:\Program Files\Carestream\Film Composer\FilmComposer.exe

C:\Program Files\Carestream\Film Composer\storescu.exe

C:\Program Files\Carestream\Film Composer\UninstallFilmComposer.exe

C:\Program Files\Carestream\Quick Burn\QuickBurn.exe

C:\Program Files\Carestream\Quick Burn\Uninstall-QuickBurn-Carestream.exe

C:\Program Files\Carestream\Report Editor\ReportEditor-uninst.exe

C:\Program Files\Carestream\Report Editor\ReportEditor.exe

C:\Program Files\Carestream\Volume Converter\VolumeConverter-uninst.exe

C:\Program Files\Carestream\Volume Converter\VolumeConverter.exe

C:\Program Files (x86)\Carestream\CSI8SuiteUninst.exe

C:\Program Files (x86)\Carestream\CSImagingClientUninst.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\config_tool.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\csds.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\CSImagingServer.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\CsisDbBackup.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\dicom_purge_tool.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\repo_check_tool.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\uninst.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\Installs\SQL\SETUP.EXE

C:\Program Files (x86)\Carestream\CS Imaging
Server\Installs\SQL\redist\DotNetFrameworks\DOTNETFX40_FULL_SETUP.EXE

C:\Program Files (x86)\Carestream\CS Imaging
Server\Installs\SQL\redist\DotNetFrameworks\DOTNETFX46-X86-X64-ALLOS.EXE

C:\Program Files (x86)\Carestream\CS Imaging
Server\Installs\SQL\redist\VisualStudioShell\VCRuntimes\VCRUNTIME140_X64.EXE

C:\Program Files (x86)\Carestream\CS Imaging
Server\Installs\SQL\x64\FIXSQLREGISTRYKEY_X64.EXE

C:\Program Files (x86)\Carestream\CS Imaging
Server\Installs\SQL\x64\FIXSQLREGISTRYKEY_X86.EXE

C:\Program Files (x86)\Carestream\CS Imaging Server\Installs\SQL\x64\LANDINGPAGE.EXE

C:\Program Files (x86)\Carestream\CS Imaging Server\Installs\SQL\x64\RSETUP.EXE

C:\Program Files (x86)\Carestream\CS Imaging Server\Installs\SQL\x64\SCENARIOENGINE.EXE

C:\Program Files (x86)\Carestream\CS Imaging Server\Installs\SQL\x64\Setup\MSMPISETUP.EXE

C:\Program Files (x86)\Carestream\CS Imaging Server\ServiceTool\service_tool.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\ServiceTool\modules\repair_tool.exe

C:\Program Files (x86)\Carestream\CS Imaging Server\ServiceTool\modules\vcredist_x86.exe

C:\Program Files (x86)\Carestream\CSDML\AcqTwain.exe

C:\Program Files (x86)\Carestream\CSDML\command_tool.exe

C:\Program Files (x86)\Carestream\CSDML\CSDMLewdog.exe

C:\Program Files (x86)\Carestream\CSDML\CSDMLite.exe

C:\Program Files (x86)\Carestream\CSDML\export_ps.exe

C:\Program Files (x86)\Carestream\CSDML\fttfx.exe

C:\Program Files (x86)\Carestream\CSDML\set_configuration.exe

C:\Program Files (x86)\Carestream\CSDML\shutdown_tool.exe

C:\Program Files (x86)\Carestream\CSDML\sqlite3.exe

C:\Program Files (x86)\Carestream\CSDML\uninstallCSDML.exe

C:\Program Files (x86)\Carestream\CSDML\IngestCodebase\codebase_cli.exe

C:\Program Files (x86)\Carestream\CSDML\IngestCodebase\ingestcodebase.exe

C:\Program Files (x86)\Carestream\Csi8WebServer-Nginx\nginx.exe

C:\Program Files (x86)\Carestream\Csi8WebServer-Nginx\NginxService.exe

C:\Program Files (x86)\Carestream\CSImaging8\AnnotationRenderer.exe

C:\Program Files (x86)\Carestream\CSImaging8\DenturesCbctToMesh.exe

C:\Program Files (x86)\Carestream\CSImaging8\QuickBurn.exe

C:\Program Files (x86)\Carestream\CSImaging8\ReconstructionHandler.exe

C:\Program Files (x86)\Carestream\CSImaging8\SetDentureProgram.exe

C:\Program Files (x86)\Carestream\CSImaging8\SetDentureProgram8x00.exe

C:\Program Files (x86)\Carestream\CSImaging8\TConvert.exe

C:\Program Files (x86)\Carestream\CSImaging8\TW.exe

C:\Program Files (x86)\Carestream\CSImaging8\TWV.exe

C:\Program Files (x86)\Carestream\CSImaging8\uilog.exe

C:\Program Files (x86)\Carestream\CSImaging8\uninst.exe

C:\Program Files (x86)\Carestream\CSImaging8\UninstallTConvert.exe

C:\Program Files (x86)\Carestream\CSImaging8\ADB2CPubClient\ADB2CPubClient.exe

C:\Program Files (x86)\Carestream\CSImaging8\FMSEditor\FMSEditor.exe
C:\Program Files (x86)\Carestream\CSImaging8\RJViewer\RJViewer.exe
C:\Program Files (x86)\Carestream\CsisMigrationTool\CsisMigrationTool.exe
C:\Program Files (x86)\Carestream\CsisMigrationTool\codebase_ingest_tool\codebase_cli.exe
C:\Program Files (x86)\Carestream\CsisMigrationTool\codebase_ingest_tool\codebase_cli_mlv.exe
C:\Program Files
(x86)\Carestream\CsisMigrationTool\codebase_ingest_tool\codebase_ingest_tool.exe
C:\Program Files (x86)\Carestream\CsisMigrationTool\jrn_ingest_tool\jrn_ingest_tool.exe
C:\Program Files (x86)\Carestream\DataCollector\ConfigDialog.exe
C:\Program Files (x86)\Carestream\DataCollector\DataCollector-uninst.exe
C:\Program Files (x86)\Carestream\DataCollector\DataCollector.exe
C:\Program Files (x86)\Carestream\DataCollector\DbMigration.exe
C:\Program Files (x86)\Carestream\DataCollector\TransactionDialog.exe
C:\Program Files (x86)\Carestream\Film Composer\FilmComposer.exe
C:\Program Files (x86)\Carestream\Film Composer\storescu.exe
C:\Program Files (x86)\Carestream\Film Composer\UninstallFilmComposer.exe
C:\Program Files (x86)\Carestream\Monitor\monitor.exe
C:\Program Files (x86)\Carestream\Monitor\UninstallMonitor.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\AcqConfigEditor.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\acquisition.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\acq_process_Win32.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\acq_process_x64.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\CSTwacker_32.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\hostwx.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\UninstallAcqTools.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\UninstallCEPHTwain.exe
C:\Program Files (x86)\Common Files\Trophy\Acquisition\UninstallPANOTwain.exe
C:\Program Files (x86)\Common Files\Trophy\CS Export\csexport.exe
C:\Program Files (x86)\Common Files\Trophy\CS Export\uninst.exe
C:\Program Files (x86)\Common Files\Trophy\CS Export\dependencies\CSCSetup.exe

C:\Program Files (x86)\Common Files\Trophy\ExtraoralFms\extraoral_fms.exe
C:\Program Files (x86)\Common Files\Trophy\ExtraoralFms\uninst.exe
C:\Program Files (x86)\Common Files\Trophy\IPET\bin\ipet.exe
C:\Program Files (x86)\Common Files\Trophy\IPET\bin\uninst.exe
C:\Program Files (x86)\Common Files\Trophy\Licensing\activate.exe
C:\Program Files (x86)\Common Files\Trophy\Licensing\activate_console.exe
C:\Program Files (x86)\Common Files\Trophy\Licensing\Uninstall activate.exe
C:\Program Files (x86)\Common Files\Trophy\Licensing\resources\elevate.exe
C:\Program Files (x86)\Common Files\Trophy\Mesh Converter\CSMeshConverter.exe
C:\Program Files (x86)\Common Files\Trophy\Mesh Converter\UninstallCSMeshConverter.exe
C:\Program Files (x86)\Common Files\Trophy\Processing\DicomConvertToFp.exe
C:\Program Files (x86)\Common Files\Trophy\Processing\DicomImg.exe
C:\Program Files (x86)\Common Files\Trophy\Processing\DicomPack.exe
C:\Program Files (x86)\Common Files\Trophy\Processing\DicomPdf.exe
C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DEclipse\EclipseUninstaller.exe
C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DEquinox\EquinoxUninstaller.exe
C:\Program Files (x86)\Common
Files\Trophy\Processing\processings2DPresets\Factory\FactoryPresetUninstaller.exe
C:\Program Files (x86)\Common Files\Trophy\RA_SDK\Appli DemoRA_SDK.exe
C:\Program Files (x86)\Common Files\Trophy\RA_SDK\Light_RA_SDKUninstaller.exe
C:\Program Files (x86)\Common
Files\Trophy\SDK\Processings2D\bin\win32\release\Processings2DSampleAdvanced.exe
C:\Program Files (x86)\Common
Files\Trophy\SDK\Processings2D\bin\win32\release\Processings2DSampleBasic.exe
C:\Program Files (x86)\Common Files\Trophy\Services\AcqNetConfService.exe
C:\Program Files (x86)\Common
Files\Trophy\Services\TaskQueue\3D_Reconstruction_SDKUninstall.exe
C:\Program Files (x86)\Common Files\Trophy\Services\TaskQueue\Sample.exe
C:\Program Files (x86)\Common Files\Trophy\Services\TaskQueue\UninstallTaskQueue.exe
C:\Program Files (x86)\Common Files\Trophy\Trace\Trace.exe
C:\Program Files (x86)\Common Files\Trophy\Trace\UninstallTrace.exe

C:\Program Files (x86)\Microsoft SQL Server\140\Shared\Sqldumper.exe
C:\Program Files (x86)\Microsoft SQL Server\140\Shared\SqIwtsn.exe
C:\Program Files (x86)\Microsoft SQL Server\140\Tools\Binn\SQLPS.exe
C:\Program Files (x86)\Microsoft SQL Server\90\Shared\sqlbrowser.exe
C:\Program Files\Microsoft SQL Server\140\COM\DISTRIB.exe
C:\Program Files\Microsoft SQL Server\140\COM\logread.exe
C:\Program Files\Microsoft SQL Server\140\COM\qrdrsvc.exe
C:\Program Files\Microsoft SQL Server\140\COM\replmerg.exe
C:\Program Files\Microsoft SQL Server\140\COM\snapshot.exe
C:\Program Files\Microsoft SQL Server\140\COM\tablediff.exe
C:\Program Files\Microsoft SQL Server\140\DTS\Binn\DTExec.exe
C:\Program Files\Microsoft SQL Server\140\DTS\Binn\dtshost.exe
C:\Program Files\Microsoft SQL Server\140\DTS\Binn\DTSWizard.exe
C:\Program Files\Microsoft SQL Server\140\DTS\Binn\dtutil.exe
C:\Program Files\Microsoft SQL Server\140\Setup Bootstrap\SQL2017\setup.exe
C:\Program Files\Microsoft SQL Server\140\Setup
Bootstrap\SQL2017\x64\FixSqlRegistryKey_x64.exe
C:\Program Files\Microsoft SQL Server\140\Setup
Bootstrap\SQL2017\x64\FixSqlRegistryKey_x86.exe
C:\Program Files\Microsoft SQL Server\140\Setup Bootstrap\SQL2017\x64\LandingPage.exe
C:\Program Files\Microsoft SQL Server\140\Setup Bootstrap\SQL2017\x64\RSetup.exe
C:\Program Files\Microsoft SQL Server\140\Setup Bootstrap\SQL2017\x64\ScenarioEngine.exe
C:\Program Files\Microsoft SQL Server\140\Setup Bootstrap\SQL2017\x64\SetupARP.exe
C:\Program Files\Microsoft SQL Server\140\Shared\Sqldumper.exe
C:\Program Files\Microsoft SQL Server\140\Shared\SqIwtsn.exe
C:\Program Files\Microsoft SQL Server\140\Shared\VS2008\1033\rdbgsetup.exe
C:\Program Files\Microsoft SQL Server\140\Tools\Binn\OSQL.EXE
C:\Program Files\Microsoft SQL Server\140\Tools\Binn\SQLdiag.exe
C:\Program Files\Microsoft SQL Server\140\Tools\Binn\SqLogShip.exe
C:\Program Files\Microsoft SQL Server\90\Shared\sqlwriter.exe
C:\Program Files\Microsoft SQL Server\Client SDK\ODBC\130\Tools\Binn\bcp.exe
C:\Program Files\Microsoft SQL Server\Client SDK\ODBC\130\Tools\Binn\SQLCMD.EXE

C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\BackupToUrl.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\DatabaseMail.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\DCEXEC.EXE
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\fdhost.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\fdlauncher.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\SQLAGENT.EXE
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\sqlceip.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\SQLIOSIM.EXE
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\sqlservr.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\StretchCodeGen.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\xpdsi.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\Xtp\VC\bin\cl.exe
 C:\Program Files\Microsoft SQL Server\MSSQL14.CSISSEVER\MSSQL\Binn\Xtp\VC\bin\link.exe

Windows Firewall

80	TCP	In/out
443	TCP	In/out
9002	TCP	In/out
9005	TCP	In/out
9999	TCP	In/out
10001-10100	TCP	In/out
50000	TCP	In/out
50920	UDP	In/out
52540	TCP	In/out
52541	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CS8100

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Common Files\Trophy\Acquisition\AcqConfigEditor.exe

C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process20_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\HostWx.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallAcqNetConfig.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallDriverCS8100.exe

C:\Program Files\Common Files\Trophy\Acquisition\Uninstall_CS8100.exe

C:\Program Files\Common Files\Trophy\Services\AcqNetConfService.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100\Application\DiagModulesConfigEditor.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100\Application\GetDISPatientPath.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100\Application\ImageQualityControl.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100\Application\ServiceTools.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100\Application\IQTools\SFR\sfr_command.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100\Application\IQTools\SFR\bin\win32\sfr2.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100\Application\IQTools\SFR\bin\win32\sfr.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100\Application\IQTools\SFR\bin\win32\sfr_command.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100\Commons\XMLTransform\xsltproc.exe

Windows Firewall

21	TCP	In/out
----	-----	--------

1001	UDP	In/out
------	-----	--------

5354	TCP	In/out
10000	TCP	In/out
10001	TCP	In/out
20000	TCP	In/out
20001	TCP	In/out
33001	TCP	In/out
33002	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CS8100SC

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Common Files\Trophy\Acquisition\AcqConfigEditor.exe

C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process20_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\HostWx.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallAcqNetConfig.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallDriverCEPH_SC.exe

C:\Program Files\Common Files\Trophy\Acquisition\Uninstall_CS8100SC.exe

C:\Program Files\Common Files\Trophy\Acquisition\Licensing\uninstall_LicensingCore.exe

C:\Program Files\Common Files\Trophy\Services\AcqNetConfService.exe

C:\Program Files\Common

Files\Trophy\ServiceTools\8100SC\Application\DiagModulesConfigEditor.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100SC\Application\GetDISPatientPath.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100SC\Application\ImageQualityControl.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100SC\Application\ServiceTools.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC\Application\IQTools\SFR\sfr_command.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC\Application\IQTools\SFR\bin\win32\sfr2.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC\Application\IQTools\SFR\bin\win32\sfr.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC\Application\IQTools\SFR\bin\win32\sfr_command.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC\Commons\XMLTransform\xsltproc.exe

Windows Firewall

21	TCP	In/out
5354	TCP	In/out
1001	UDP	In/out
10000	TCP	In/out
10001	TCP	In/out
20000	TCP	In/out
20001	TCP	In/out
33001	TCP	In/out
33002	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CS81003D

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Common Files\Trophy\Acquisition\AcqConfigEditor.exe

C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process20_x64.exe
 C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe
 C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe
 C:\Program Files\Common Files\Trophy\Acquisition\HostWx.exe
 C:\Program Files\Common Files\Trophy\Acquisition\UninstallAcqNetConfig.exe
 C:\Program Files\Common Files\Trophy\Acquisition\UninstallDriverCS81003D.exe
 C:\Program Files\Common Files\Trophy\Acquisition\Uninstall_CS81003D.exe
 C:\Program Files\Common Files\Trophy\Services\AcqNetConfService.exe
 C:\Program Files\Common Files\Trophy\Services\TaskQueue\Sample.exe
 C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePanel.exe
 C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePluginProcess.exe
 C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueueServer.exe
 C:\Program Files\Common Files\Trophy\Services\TaskQueue\Tools\CudaInfo.exe
 C:\Program Files\Common
 Files\Trophy\ServiceTools\8100_3D\Application\DiagModulesConfigEditor.exe
 C:\Program Files\Common Files\Trophy\ServiceTools\8100_3D\Application\GetDISPatientPath.exe
 C:\Program Files\Common
 Files\Trophy\ServiceTools\8100_3D\Application\ImageQualityControl.exe
 C:\Program Files\Common Files\Trophy\ServiceTools\8100_3D\Application\ServiceTools.exe
 C:\Program Files\Common
 Files\Trophy\ServiceTools\8100_3D\Application\IQTools\SFR\sfr_command.exe
 C:\Program Files\Common
 Files\Trophy\ServiceTools\8100_3D\Application\IQTools\SFR\bin\win32\sfr2.exe
 C:\Program Files\Common
 Files\Trophy\ServiceTools\8100_3D\Application\IQTools\SFR\bin\win32\sfr.exe
 C:\Program Files\Common
 Files\Trophy\ServiceTools\8100_3D\Application\IQTools\SFR\bin\win32\sfr_command.exe
 C:\Program Files\Common
 Files\Trophy\ServiceTools\8100_3D\Commons\XMLTransform\xsltproc.exe

Windows Firewall

5354	TCP	In/out
1001	UDP	In/out
10000	TCP	In/out
10001	TCP	In/out
20000	TCP	In/out
20001	TCP	In/out
33001	TCP	In/out
33002	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CS8100SC3D

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Common Files\Trophy\Acquisition\AcqConfigEditor.exe

C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process20_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\HostWx.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallAcqNetConfig.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallDriver3D_SC.exe

C:\Program Files\Common Files\Trophy\Acquisition\Uninstall_CS8100SC3D.exe

C:\Program Files\Common Files\Trophy\Services\AcqNetConfService.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Sample.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePanel.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePluginProcess.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueueServer.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Tools\CudaInfo.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC3D\Application\DiagModulesConfigEditor.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100SC3D\Application\GetDISPatientPath.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC3D\Application\ImageQualityControl.exe

C:\Program Files\Common Files\Trophy\ServiceTools\8100SC3D\Application\ServiceTools.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC3D\Application\IQTools\SFR\sfr_command.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC3D\Application\IQTools\SFR\bin\win32\sfr2.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC3D\Application\IQTools\SFR\bin\win32\sfr.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC3D\Application\IQTools\SFR\bin\win32\sfr_command.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\8100SC3D\Commons\XMLTransform\xsltproc.exe

Windows Firewall

21	TCP	In/out
5354	TCP	In/out
1001	UDP	In/out
10000	TCP	In/out
10001	TCP	In/out
20000	TCP	In/out
20001	TCP	In/out
33001	TCP	In/out
33002	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CS8200v1

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Common Files\Trophy\Acquisition\AcqConfigEditor.exe

C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process20_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\hostwx.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallAcqNetConfig.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallDriverSyracuse.exe

C:\Program Files\Common Files\Trophy\Acquisition\Uninstall_CS82003D.exe

C:\Program Files\Common
Files\Trophy\Acquisition\plugins\AcqPluginIOScan\Uninstall_AcqPluginIOScan.exe

C:\Program Files\Common Files\Trophy\Services\AcqNetConfService.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Sample.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePanel.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePluginProcess.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueueServer.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Tools\CudaInfo.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\82003D\Application\DiagModulesConfigEditor.exe

C:\Program Files\Common Files\Trophy\ServiceTools\82003D\Application\ImageQualityControl.exe

C:\Program Files\Common Files\Trophy\ServiceTools\82003D\Application\ServiceTools.exe

C:\Program Files\Common Files\Trophy\ServiceTools\82003D\Application\TrackBalls_UnitTests.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\82003D\Application\IQTools\SFR\sfr_command.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\82003D\Application\IQTools\SFR\bin\win32\sfr2.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\82003D\Application\IQTools\SFR\bin\win32\sfr.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\82003D\Application\IQTools\SFR\bin\win32\sfr_command.exe

C:\Program Files\Common
Files\Trophy\ServiceTools\82003D\Commons\XMLTransform\xsltproc.exe

Windows Firewall

21	TCP	In/out
5354	TCP	In/out
1001	UDP	In/out
10000	TCP	In/out
10001	TCP	In/out
20000	TCP	In/out
20001	TCP	In/out
33001	TCP	In/out
33002	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CS8200v2

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Common Files\Trophy\Acquisition\AcqConfigEditor.exe

C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process20_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\hostwx.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallAcqNetConfig.exe

C:\Program Files\Common Files\Trophy\Acquisition\UninstallAcquisitionGui_Syracuse.exe

C:\Program Files\Common Files\Trophy\Acquisition\Uninstall_CS82003D.exe

C:\Program Files\Common
Files\Trophy\Acquisition\DeviceStateAgent\UninstallDeviceStateAgent_Radiology.exe

C:\Program Files\Common
Files\Trophy\Acquisition\DeviceStateAgent\Radiology\Application\DeviceStateAgent.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServer\GuiServer.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServer\HkcuRegkeyAccess.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServer\UninstallGuiServerSyracuseV2.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServerService\GUIServerService.exe

C:\Program Files\Common
Files\Trophy\Acquisition\GUIServerService\Uninstall_SetupGUIServerSyracuseV2Service.exe

C:\Program Files\Common Files\Trophy\Acquisition\Licensing\uninstall_LicensingCore.exe

C:\Program Files\Common Files\Trophy\Acquisition\NagoyaWebView\Nagoya_WebView.exe

C:\Program Files\Common
Files\Trophy\Acquisition\NagoyaWebView\Uninstall_NagoyaWebView.exe

C:\Program Files\Common Files\Trophy\Acquisition\nginx\nginx.exe

C:\Program Files\Common Files\Trophy\Acquisition\nginx\NginxParserConf.exe

C:\Program Files\Common Files\Trophy\Acquisition\NGinxService\NGinxService.exe

C:\Program Files\Common Files\Trophy\Acquisition\NGinxService\Uninstall_SetupNGinxService.exe

C:\Program Files\Common
Files\Trophy\Acquisition\PluginBackEnds\AcqSyracuseBackEnd\Configure.exe

C:\Program Files\Common
Files\Trophy\Acquisition\PluginBackEnds\AcqSyracuseBackEnd\UninstallDriverSyracuseBackEnd.exe

C:\Program Files\Common
Files\Trophy\Acquisition\Plugins\AcqPluginIOScan\Uninstall_AcqPluginIOScan.exe

C:\Program Files\Common Files\Trophy\Acquisition\Plugins\AcqSyracuse\Configure.exe

C:\Program Files\Common
Files\Trophy\Acquisition\Plugins\AcqSyracuse\Uninstall_DriverSyracuse.exe

C:\Program Files\Common Files\Trophy\Acquisition\WSSurveyServer\UninstallWSSurveyServer.exe

C:\Program Files\Common Files\Trophy\Acquisition\WSSurveyServer\WSSurveyServer.exe

C:\Program Files\Common
Files\Trophy\Acquisition\WSSurveyService\Uninstall_SetupWSSurveyService_Syracuse.exe

C:\Program Files\Common Files\Trophy\Acquisition\WSSurveyService\WSSurveyReport.exe

C:\Program Files\Common Files\Trophy\Acquisition\WSSurveyService\WSSurveyService.exe

C:\Program Files\Common Files\Trophy\AEC\aec.exe

C:\Program Files\Common Files\Trophy\AEC\CreatesListProjectionMatricesFromDCM.exe

C:\Program Files\Common Files\Trophy\AEC\EXE_CreateListProjectionMatrices.exe

C:\Program Files\Common Files\Trophy\AEC\vcredist2013_x64.exe

C:\Program Files\Common Files\Trophy\AEC\vcredist_x64.exe

C:\Program Files\Common Files\Trophy\AEC\VC_redist.x64.exe

C:\Program Files\Common Files\Trophy\AEC\SFR\sfr_command.exe

C:\Program Files\Common Files\Trophy\AEC\SFR\bin\win32\sfr2.exe

C:\Program Files\Common Files\Trophy\AEC\SFR\bin\win32\sfr.exe

C:\Program Files\Common Files\Trophy\AEC\SFR\bin\win32\sfr_command.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\AutoStart.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWUninstall.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\Etk.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\Install.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\PP_Update.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\Uninstall_CS8200EParts.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\h2\service\wrapper.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\jabswitch.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\java-rmi.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\java.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\javacpl.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\javaw.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\javaws.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\jjs.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\jp2launcher.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\keytool.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\kinit.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\klist.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\ktab.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\orbd.exe

C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\pack200.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\policytool.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\rmid.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\rmiregistry.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\servertool.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\ssvagent.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\tnameserv.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\32\bin\unpack200.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\jabswitch.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\java-rmi.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\java.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\javacpl.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\javaw.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\javaws.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\jjs.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\jp2launcher.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\keytool.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\kinit.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\klist.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\ktab.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\orbd.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\pack200.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\policytool.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\rmid.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\rmiregistry.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\servertool.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\ssvagent.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\tnameserv.exe
C:\Program Files\Common Files\Trophy\CS8200 eParts\DWJava\Java8u131\64\bin\unpack200.exe
C:\Program Files\Common Files\Trophy\EquipmentTools\UninstallEquipmentTools_CS8200.exe
C:\Program Files\Common Files\Trophy\EquipmentTools\CS8200\Application\EquipmentTools.exe

C:\Program Files\Common
Files\Trophy\EquipmentTools\CS8200\Application\IQTools\SFR\sfr_command.exe

C:\Program Files\Common
Files\Trophy\EquipmentTools\CS8200\Application\IQTools\SFR\bin\win32\sfr2.exe

C:\Program Files\Common
Files\Trophy\EquipmentTools\CS8200\Application\IQTools\SFR\bin\win32\sfr.exe

C:\Program Files\Common
Files\Trophy\EquipmentTools\CS8200\Application\IQTools\SFR\bin\win32\sfr_command.exe

C:\Program Files\Common
Files\Trophy\EquipmentTools\CS8200\Commons\XMLTransform\xsltproc.exe

C:\Program Files\Common Files\Trophy\Services\AcqNetConfService.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Sample.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePanel.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePluginProcess.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueueServer.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Tools\CudaInfo.exe

Windows Firewall

21	TCP	In/out
5354	TCP	In/out
1001	UDP	In/out
10000	TCP	In/out
10001	TCP	In/out
20000	TCP	In/out
20001	TCP	In/out
33001	TCP	In/out
33002	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

CS9600

Windows Defender

Folders:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream

Executables:

C:\Program Files\Common Files\Trophy\Acquisition\acquisition.exe

C:\Program Files\Common Files\Trophy\Acquisition\AcquisitionSampleAdvanced.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process20_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_Win32.exe

C:\Program Files\Common Files\Trophy\Acquisition\acq_process_x64.exe

C:\Program Files\Common Files\Trophy\Acquisition\Uninstall_CS9600Proxy.exe

C:\Program Files\Common
Files\Trophy\Acquisition\DeviceStateAgent\UninstallDeviceStateAgent_Dental.exe

C:\Program Files\Common
Files\Trophy\Acquisition\DeviceStateAgent\Dental\Application\DeviceStateAgent.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\Uninstall_SetupGStreamer.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gdbus.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gdk-pixbuf-csource.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gdk-pixbuf-query-loaders.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gio-querymodules.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\glib-compile-schemas.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gsettings.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gst-device-monitor-1.0.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gst-discoverer-1.0.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gst-inspect-1.0.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gst-launch-1.0.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gst-play-1.0.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\gst-typefind-1.0.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\json-glib-format.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\bin\json-glib-validate.exe

C:\Program Files\Common Files\Trophy\Acquisition\GStreamer\libexec\gststreamer-1.0\gst-plugin-scanner.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServer\GuiServer.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServer\HkcuRegkeyAccess.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServer\UninstallGuiServerSyracuseV1.exe

C:\Program Files\Common Files\Trophy\Acquisition\GUIServerService\GUIServerService.exe

C:\Program Files\Common
Files\Trophy\Acquisition\GUIServerService\Uninstall_SetupGUIServerSyracuseV1Service.exe

C:\Program Files\Common Files\Trophy\Acquisition\NagoyaWebView\Nagoya_WebView.exe

C:\Program Files\Common
Files\Trophy\Acquisition\NagoyaWebView\Uninstall_NagoyaWebView.exe

C:\Program Files\Common Files\Trophy\Acquisition\nginx\nginx.exe

C:\Program Files\Common Files\Trophy\Acquisition\nginx\NginxParserConf.exe

C:\Program Files\Common Files\Trophy\Acquisition\NGinxService\NGinxService.exe

C:\Program Files\Common Files\Trophy\Acquisition\NGinxService\Uninstall_SetupNGinxService.exe

C:\Program Files\Common
Files\Trophy\Acquisition\plugins\AcqNagoyaProxy\Uninstall_DriverCS9600Proxy.exe

C:\Program Files\Common
Files\Trophy\Acquisition\plugins\AcqPluginIOScan\Uninstall_AcqPluginIOScan.exe

C:\Program Files\Common
Files\Trophy\Acquisition\WorkStationDiscovery\UninstallWorkStationDiscovery.exe

C:\Program Files\Common
Files\Trophy\Acquisition\WorkStationDiscovery\WorkStationDiscovery_Client.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Sample.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePanel.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueuePluginProcess.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\TaskQueueServer.exe

C:\Program Files\Common Files\Trophy\Services\TaskQueue\Tools\CudaInfo.exe

Windows Firewall

80	TCP	In/out
443	TCP	In/out
3389	TCP	In/out
5354	TCP	In/out

Windows Security

Add rights to the current user, "Read and Execute" on the folders below:

C:\Program Files\Carestream

C:\Program Files\Common Files\Trophy

C:\Program Files (x86)\Carestream

C:\Program Files (x86)\Common Files\Carestream

C:\Program Files (x86)\Common Files\Trophy

Add rights to the current user, "Full Control" on the folders below:

C:\ProgramData\TW

C:\Users\Public\Documents\Carestream